NASO BILLIARY DRAINAGE CATHETER



Key features:

- Pigtail shape for better retention.
- Each set includes nasal transfer tube and a connector to facilitate irrigation and drainage
- Available as Straight, Straight with Pigtail and Angled with Pigtail variants

Specifications

Size

Diameter 5, 6, 7, 8, 8.5 & 10 Fr

Length 250 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acmedicaldevices.com, acmedicaldevices.com