NASO JEJUNAL FEEDING TUBE

Used to provide short duration enteral access for delivery of nutrition or medication to the small bowel. Supplied sterile in a peel open Tyvek pouch. Designed for single use only.



Key features:

- Smooth surface for better patient comfort
- Accepts 0.035 inch guidewire
- Supplied with or without guidewire

Specifications

Size

Diameter 08, 10, 12, 14, 16, 18 Fr

Length 120, 155 and 200 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com