NEEDLE ELECTRODE - BIPOLAR

Needle Electrode is used to make Plasma incisions during TURP Surgery. Product supplied Sterile and single use only.



Key features:

- Faster firing, longer lasting
- High quality insulation.
- .Compatible with all standard resectoscopes.
- Available in 30 and 45 degrees.

Classification:

Compatible with

- STORZ 1 Stem Instrument
- STORZ 2 Stem Instrument
- Olympus 1 Stem Instrument
- Wolf 1 Stem Instrument

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

Other properties:

Available with / without bipolar cable.

Specifications

Size

Diameter 24
Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com