NOTTINGHAM ONE STEP DILATOR

Used to dilate the intramural tunnel and ureter in one step. May also be used to inject contrast medium to determine sub mucosal tunneling or anesthetic jelly to aid in dislodging ureteral stones. Designed for single use only.



Key features:

- Smooth surfaces for ease of introduction
- Uniform taper to reduce trauma.
- Used with 0.038" guide wire.

Specifications

Size

Diameter 8/4, 9/5, 10/5, 10/6, 10/7, 12/6, 12/7, 14/6, 14/7, 16/6, 16/7, 16/8 and 18/6 Fr

Length 70 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com