
PCN SET - MALECOT

Used for Percutaneous placement of a malecot catheter in the renal pelvis for nephrostomy drainage. Designed for single use only.



Key features:

- Highly radiopaque for better visulization.
- The Malecot flower design enhances drainage and promotes catheter retention.
- Round tip malecot catheters reduces the trauma to the tissues.
- Highly suited for drainage of thick viscous fluids.
- Malecot Catheter can be supplied with / without Hydrophilic coating

Classification:

PCN Malecot Set consists of:

- Malecot Catheter
- Set of Fascial Dilators
- Initial Puncture Needle
- "J" Tip Guidewire

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

- Urine Bag Connector
- Scalpel

Available in Tray packaging as well as pouch packaging.

Specifications

Size

Diameter 6, 8, 10, 12, 14, 16, 18, 20,22 Fr

Length 20 - 35 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com