STONE BASKET - PERK BASKET (PCNL)

Perk Baskets are used for kidney stone removal, under direct vision. Intended for single use only.



Key features:

- Detachable basket
- Nitinol wire Tipless basket

Specifications

Size

Diameter 10, 12 Fr

Length 37 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com