PTBD CATHETER

Used for Percutaneous transhepatic biliary drainage (PTBD). Designed for single use only.



Key features:

- Pigtail Shape for better retention
- Smooth surface for trauma
- Available as straight, angled, angled with pigtail
- Available as catheter and set
- Set includes PTBD catheter, Dilators, needle and guidewire

Specifications

Size

Diameter 6,7,8,8.5,10,12 and 14 Fr

Length 35, 40, 45 and 60 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com