
S Shaped Dilator Set

Used for dilation of male urethral strictures. It supports the natural curvature of a male urethra. Designed for single use only.



Key features:

- Smooth surfaces for ease of introduction
- Uniform taper to reduce trauma
- Dilators provided with hydrophilic coating
- Radioaque for better fluoroscopic visualisation
- Used with 0.038" guide wire.

Specifications

Size

Diameter 8-20 Fr

Length 37 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com