SCLEROTHERAPY / INJECTION NEEDLE

Used for endoscopic injection into gastrointestinal mucosa. Supplied sterile in a peel open Tyvek pouch. For single use only.



Key features:

- Efficient, smooth and easy to use handle
- Exceptional needle actuation even in tortuous pathways
- Low coefficient of friction for smooth and consistent flow
- Kink resistant tubing

Specifications

Size

Diameter 2.7 mm

Length 1600,1700, 1750, 1800, 2000, 2200, 2250 and 2300 mm

Needle Size 21, 23 and 25 G

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com