
SPHINCTEROTOME / PAPILLOTOME

Used for Endoscopic Cannulation of Biliary Ductal System and Sphincterotomy.



Key features:

- Triple Lumen with pre curved tapered Tip
- Cutting Wire Length: 20, 30 mm
- Tip design: 5.5 Fr
- Available with / without insulation
- Efficient, smooth and easy to use handle
- Recommended Guidewire 0.035 inch.

Specifications

Size

Sheath Diameter 6.8 Fr

Length 1800 and 2000

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

mm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com