STONE BASKET -STONE HOLD

Used to capture and remove calculi and other foreign bodies in the urinary tract, and to minimize stone migration during intracorporeal lithotripsy. Intended for single use only.



Key features:

- Prevents stone migration during lithotripsy
- Eliminates the need for flexible scopes during procedure
- Ergonomic design

Specifications

Size

Diameter 1.9, 2.2, 2.4 and 3 Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com