STONE RETRIEVAL BASKET NITINOL

Intended to retrieve Foreign Bodies and C.B.D. stones endoscopically. Designed for single use only.



Key features:

- Engineered to facilitate pushability
- Unique handle for ease of use
- Atraumatic tip for ease of procedure
- Compatible with 2.8 mm instrument channel
- Hexagonal, Trapezoidal and Balloon shaped

Specifications

Size

Diameter 15, 20, 25 30 and 35 mm

Length 2300 mm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com