STONE BASKET - SUPER CLASP NITINOL STONE RETRIEVER

Designed for the manipulation and retrieval of stones from the Bladder, Ureter and Renal Pelvis. This device is supplied sterile in a peel-open pouch. Intended for single use only.



Key features:

- Prevents retropulsion of stones up the ureter and into the kidney
- Distinctive and slim design
- Sturdy and kink resistant
- Choice of handles provided

Specifications

Size

Diameter 1.9, 2.2, 2.4 and 3

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

Fr

Length 70,90,110,120,130 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com