SUPRAPUBIC MALECOT CATHETER (SPC)

Used to provide nephrostomy drainage by percutaneous placement of a malecot catheter in the renal pelvis.



Key features:

- Highly radiopaque for better visualization
- The Malecot flower design enhanced drainage & promote catheter retention
- Round tip malecot catheters reduces the trauma to the tissue

Classification:

Malecot Catheter supplied with a metal trocar.

Specifications

Size

Diameter 6, 8, 10, 12, 14, 16, 18, 20,22 Fr

Length 20 - 40 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com