STONE BASKETS - SUPER CLENCH NITINOL RETRIEVER

Designed for the manipulation and retrieval of stones from the Bladder, Ureter and Renal Pelvis. This device is supplied sterile in a peel-open pouch. Intended for one-time use.



Key features:

- Multiwire design for tight clenching of fragments
- Partial closure forms a web of wires preventing retropulsion of stones
- Sturdy and kink resistant
- Unique design

Specifications

Size

Diameter 1.9,2.2,2.4 and 3 Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

Length	70,90,110,120 and 130 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com