
# TRI PRONG GRASPER

Used to facilitate endoscopic manipulation/ removal of stones and retrival of foreign objects. Designed for single use only.



## Key features:

- Designed for minimum trauma.
- Maximum grasping effciency.

## Classification:

For URS and Cystoscope

- Rigid
- Flexible
- PCNL

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

# Other properties:

Rigid

Size : 1.9, 2.5, 3, 4, 5 and 6

Length : 40 and 60 cms

Flexible

Size : 1.9, 2.5, 3, 4, 5 and 6

Length : 40 and 60 cms

PCNL

Size : 7, 9 and 12 Fr

Length : 40 cms

## Specifications

Size

**Diameter** 1.9 - 12 Fr

**Length** 40 - 60 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com