URETERAL ACCESS SHEATH

Used for ureteral dilation and to establish a conduit during endoscopic urological procedures facilitating the passage of flexible scopes into the urinary tract. Designed for single use only



Key features:

- Atraumatic tip for easy access
- Kink resistant coil for smooth navigation through a torturous anatomy
- Irradium-Platinum marker for highest radiopacity
- Tapered dilator for easy intramural access
- Can be supplied with hydrophilic coating

Specifications

Size

Diameter 9.5/11.5, 10.5/12.5, 12/14 Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

Length 28, 35, 38, 45 and 55 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com