
URETERAL DILATOR SET

Used for dilation of the ureter prior to uretero renoscopy and also for sequential (step by step) dilation of the intramural tunnel and ureter and/or stone manipulation. Designed for single use only.



Key features:

- Smooth surfaces for ease of introduction
- Uniform taper to reduce trauma
- Can be supplied with / without hydrophilic coating

Specifications

Size

Diameter 6 - 16, 6- 22 Fr

Length 70 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com