# URETHRAL DILATOR SET

Used for dilation of the urethral strictures and to establish and maintain an access in to the bladder via urethra for cystoscopy. Designed for single use only.



## Key features:

- Smooth surfaces for ease of introduction
- Special atraumatic design
- Radiopaque for clear fluoroscopic visualisation

## Specifications

Size

**Diameter** 8, 10, 12, 14, 16, 18, 20, 22 and 24 Fr

**Length** 40 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com